

Drug-Free Antibacterial Hybrid Biopolymers for Medical Applications

S. aureus Colonized Human Skin Equivalent - In Vitro Bio-evaluation Tool for **Antibacterial Polymers: ESR14**

Ayesha Idrees^{1, 2}, Richard Viebahn², Gianluca Ciardelli¹, Valeria Chiono¹, Jochen Salber²

¹Department of Mechanical and Aerospace Engineering (DIMEAS), Politecnico di Torino, Italy

² UMC Knappschaftskrankenhaus GmbH, Clinic of Surgery, Hospital of the Ruhr-University, Bochum, Germany


A bacterial colonized human skin equivalent (c-HSE)

The aim of this study was the development of a human skin model and human skin wound infection model for the bio-evaluation of antimicrobial biomaterials intended for wound healing purposes. The three-dimensional in vitro models will represent advanced and complex systems to perform more reliable preclinical studies. These models will be employed for *in vitro* screening of both antibacterial activity as well as cytocompatibility of new biomaterials and wound dressings to optimize their in vivo performances. In the study, the 3D systems were developed and their structure was characterized. The antibacterial activity and cytotoxicity of a model wound dressing releasing Ag⁺ was analyzed in the models.






3D Dermal Fibroblast Model



Optimizing the dermal part of human skin: Z-stacked imaging revealed the filopodia like morphology and a uniform distribution of human fibroblasts at different planes inside a Col-I matrix. Fluorescent microscopic images show cell nuclei stained with DAPI and cytoskeletal F-actin stained with Phalloidin. Scale bar=100 µm

Histological Analysis of the HSE

Immunohistochemistry of the 3D Human Skin Model




S. aureus-colonized HSE

H&E stained cross

TEM image:

Ultrastructure Analysis


Inoculated bacteria adhere to the




section of *in vitro* HSE model:

Epidermal layer (1) Dermal layer (2)

The HSE epidermis has a characteristic structure: Stratum corneum, granulosum, spinosum and basale.



Collagen-I fibres (Col I) Epidermal-**Dermal Junction** (DEJ). **DEJ** presents lamina lucida (LL) and lamina densa (LD).




dermal surface, colonize, and replicate to make large structures of biofilm.

Big arrows: Bacteria located within a biofilm matrix inside the dermis.

Small arrows: Bacteria surrounding keratinocytes in epidermis.

Cytocompatibility Analysis



(a) (b) reus in (%) 9.8 100 Colonize Log) 9.4 au 80 60 8.6 educti Broth ; `≦ 8.2 rcent R Liquid 8.7 De Skin 7.4 20 7.4 Log Ag+ MBC Ag+ wound Ag+ MBC Ag+ Untreated dressing

Antibacterial Analysis

Decreasing concentration of Ag+

Ag+ wound Wound Untreated dressing dressing control control (Without Ag+)

wound dressing

Cell viability measuring in the 3D system. The 3D skin model was exposed to a range of silver ion concentrations (Ag⁺) for a period of 24 hours. A commercially available Ag⁺ releasing wound dressing served the purpose of a model material and was tested in a 3D system along with its control material (without Ag⁺).

The graph demonstrates the treatment of infected skin equivalents with a commercially available Ag⁺ releasing wound dressing. Skin equivalents were infected with S.aureus and thereafter, Ag⁺ releasing wound dressing or Ag⁺ in PBS was applied onto the skin equivalents.

Conclusion: Development of colonized human skin equivalent (c-HSE); Risk assessment platform for cytocompatibility evaluation; Efficacy assessment of antibacterial materials; Comparison of 2D vs. 3D systems; Understanding "Host-Pathogen Interaction"; Development of complex skin models.



HyMedPoly has received funding from the European Union's Horizon 2020 research and innovation programme under the Marie Sklodowska-Curie grant agreement number 643050.

www.hymedpoly.eu

ESR14 Email: ayesha.idrees@polito.it; ayeshaidrees19@gmail.com Supervisors' Emails: gianluca.ciardelli@polito.it; valeria.chiono@polito.it; jochen.salber@kk-bochum.de; jochen.salber@hotmail.com Affiliation: Department of Mechanical and Aerospace Engineering (DIMEAS), Politecnico di Torino, Italy; UMC Knappschaftskrankenhaus GmbH, Clinic of Surgery, Hospital of the Ruhr-University, Bochum, Germany